


REVIVE MOTION

Complete App | All 10 Screens | MIT App Inventor





3 benefits first


Show what the app does before asking for anything. Builds motivation to complete setup.

Injury selector

Determines exercise library shown. Saves injury_type to TinyDB. 4 large tap targets.

Friendly CTA

'Let's Go!' not 'Next'. Green gradient signals positive. 5 steps total, only 2 shown.



Today's exercise

Pulled from PT prescription in CloudDB. One-tap video link. Patient knows exactly what to do.

PT update flag


If PT pushed changes, patient sees it before session. Not buried in settings.

Start button

Large, impossible to miss. App does not jump to live view — patient reviews first.

Post-session summary

After session ends: improvement vs last time. The most motivating moment in the app.



3-state feedback

Green / Amber / Red. 'Almost there!' for 88-99% of target. 46pt font readable across room.

Rep counter

Large pill shape, front and center after feedback. Counts automatically when form is correct.

Live angle chart


Real-time joint angle waveform. Gold dot = current. Green goal line at 90d.

EMG waveform

Orange zones show when correct muscle is firing. Confirms the right muscle is working.

Exercise card

PT-prescribed exercise + pain score + video demo link. All in one card.



Scale 1-10


Color coded: green (1-3), orange (4-6), red (7-10). Quick and familiar.

Clinically critical

If 7+: session blocked, caregiver CloudDB alert fired. No other student app does this.

Pain trend chart

7-day bar chart visible on the check-in screen. Caregiver sees same on their dashboard.




XP level bar
Avatar speech bubble
Fills toward MVP. Level name changes with each character peak score. Message changes based on performance.

Quest front and center
Daily quest with progress bar above the rep button. Main goal, always visible.

Tap to Complete Rep
Label tells the kid exactly what to do. +10 for correct, +20 for perfect form.

4 achievement badges
Locked/unlocked. Milestone-based: 1st rep, 3-day streak, 7-day streak, ROM goal.



ROM line chart

12 sessions. Missed sessions = red dot at baseline. Green dashed goal line at 90d.

Compliance heatmap


4-week calendar. Dark green = full session. Fading to empty = missed or partial.

Weekly reps bars

Color-coded: green=10/10, orange=partial, red=missed. Week at a glance.

Pain trend line

Falling line = recovery working. Rising line = flag for caregiver attention.



Status first
"NEEDS ATTENTION" or "ON TRACK" in large text. "Send Reminder" lives in the alert card. No navigation needed to act.

Dual-axis chart
ROM rising (blue), pain falling (dashed red). Inverse trend shows real recovery.

Before/Now/Goal
Single chart tells the whole story. Baseline to current to goal in one view.

Session calendar
Month heatmap. Missed day patterns visible instantly without counting.




4 donut rings
Compliance, ROM, weekly gain, effort. Progress rings are faster to read than numbers.

Recovery radar
5-axis. Dennis vs average Week-12 patient. Shows where he leads and where to focus.

AI analysis
3-line summary. Language toggle switches to Spanish. Generated from real session data.

Exercise prescription
Sets/ reps/hold/target as visual pills. Video link. 'Why' explanation collapsible.

Physician report
AI-written summary. Auto-sends every Sunday. Doctor reads in 30 seconds.




PT sees patient data
Compliance, ROM, and streak shown before making changes. Data-informed decisions.

Exercise dropdown
PT selects from exercise library. Sets reps, hold time, sets, and target angle.

4 metric pills
Same visual format patient sees on briefing screen. What you set is what they see.

Push to patient
One tap. CloudDB saves instantly. Patient sees PT notes on their next briefing screen.



Injury filter
Shows only exercises relevant to the injury type set during onboarding.

Prescribed badge
PT-prescribed exercise highlighted in green. Patient always knows what's required.

Video link
Every exercise has a demo video link. Tap opens in WebViewer component.

Stage conditions
Shows when each exercise becomes available based on current ROM.